Gamma Ray Spectroscopy

Abstract: A NaI(Tl) scintillation spectrometer was used to analyze the spectra of several radionuclides. We first took spectra for ¹³⁷Cs, ⁶⁰Co and ²²Na to calibrate our setup and learn about gamma decay in these isotopes. We then used gamma spectroscopy to identify the constituents of two unknown radioactive sources. Finally, we attempted to determine the practicality of using gamma spectroscopy to detect copper.

Introduction

Gamma rays are highly penetrating photons produced by positron-electron annihilation or by the decay of a radioactive nucleus. Gamma decay tends to occur following an alpha or beta decay to bring the nucleus down to ground state. These photons were first detected by Antoine Henri Becquerel in the late 1800's. Becquerel deduced that gamma rays interact with matter like light using photographic film to detect the radiation. The detection of gamma rays using scintillation spectroscopy rests on this discovery that gamma rays are photons. Gamma rays interact with matter by three processes- the photoelectric effect, Compton scattering, and pair production.

Each of these processes transfers energy to electrons, which then lose energy by ionizing atoms as they travel through matter. In the case of the photoelectric effect, all of the energy of the gamma photon is transferred to an electron. Since the ionization energy of the atoms in the scintillator is small compared to that of a gamma ray, the electron is considered to have an energy equal to that of the incident photon. This process is dominant at gamma energies lower than 511 KeV. For a nuclear charge Z, the photoelectric cross section for a K shell electron is:

$$\sigma_{photo} \approx \sigma_T Z^5 \left[\frac{h v}{m c^2} \right]^{\frac{7}{2}}$$
 (1)

Compton scattering is a demonstration of the particle-like behavior of light. During this process, a photon and an electron collide elastically, producing a scattered photon of lower energy and an electron with the energy lost by the photon. The energy of the new photon is dependent on the angle of deflection. The energies of the scattered photon (E_{ν}) and electron $(E_{\rm e})$ for an incident photon of energy E_{ν} are:

$$E'_{\gamma} = \frac{E_{\gamma}}{1 + \left(E_{\gamma} / mc^2\right)(1 - \cos\theta)} \qquad E_{e} = E'_{\gamma} - E_{\gamma}$$
(3)

The Klein-Nishina formula gives the differential cross section for Compton scattering:
$$\frac{d\sigma}{d\Omega} = \frac{1}{2} r_{\epsilon}^2 (P(E_{\gamma},\theta) - P(E_{\gamma},\theta)^2 \sin^2(\theta) + P(E_{\gamma},\theta)^3) \tag{4}$$

where $P(E_{\gamma},\theta)$ is the ratio of the final photon energy to the initial photon energy, r_e is the classical electron radius and $d\Omega$ is the solid angle. Compton scattering is dominant for middle-energy gamma rays.

At energies higher than the rest mass of two electrons ($2m_ec^2=1.022$ MeV), pair production is the principal mechanism for gamma interactions. During this process, a photon turns into an electron and a positron. This must occur near a nucleus in order to conserve momentum and energy.

All three of these processes result in the transfer of energy from a photon to an electron. As the electron travels through matter, it ionizes atoms through a series of

inelastic collisions. This ionization is measurable by a number of techniques making it useful for detecting particles.

Experimental Technique

As described above, gamma rays deposit their energy in matter such that a number of ionized atoms are created. In a scintillator, this ionization energy is converted into visible light. This property was first employed in 1909 by Ernest Rutherford who used a ZnS screen to detect scattered alpha particles. The visible light produced in the scintillating crystal is linearly proportional to the energy input of the crystal. In this experiment, we used sodium-iodide, an inorganic crystal, doped with thallium, which has an efficiency of 12% and produces 4×10^4 photons per incident MeV and has a time constant of ~200 ns. The NaI(Tl) crystal is enclosed in a sealed can and optically connected to a photomultiplier.

We used a standard photomultiplier tube (PMT) from RCA. The PMT proportionally converts the visible light from the scintillation crystal into a small electrical pulse. The first part of a PMT is a photocathode, which converts light into electrons by the photoelectric effect. Next comes a series of dynodes made of materials with good secondary electron emission, between which a high voltage is applied. This creates a potential difference of 50 to 100V between dynode pairs. The HV source used in this experiment is an Ortec 478 HV supply at positive 1.1 kV. The electrons emitted from the photocathode are accelerated to the first dynode, which then emits a proportional number of electrons to the next dynode. This signal continues to amplify along the chain of dynodes, adding up to a total amplification of $\sim 10^6$. This pulse is collected at the anode, which is connected to a preamplifier. The photomultiplier sends out a negative pulse of charge with a width on the order of 6µs and amplitude of about 200 mV.

The preamplifier converts the charge on the anode to a positive voltage signal with width of about 140µs and sends it to a pulse shaping amplifier. An ORTEC 485 amplifier was used in this experiment. The pulse shaping amplifier removes the effect of the scintillator time constant on the width of the pulse. The pulse width is now controlled by the shaping time of the amplifier. Typical pulse widths were on the order of 7µs with a height of around 25 V. The now amplified voltage signal is digitized using a Pocket Multi Channel Analyzer analog-to-digital converter. The MCA is connected to a computer in which a histogram of voltage pulse heights is plotted. This plot gives the energy spectrum of a gamma ray source, once channels are calibrated to energies. Cesium-137 was used to calibrate pulse height to photon energy (see appendix for details).

The ⁶⁰Co Spectrum

The cobalt-60 gamma spectrum shows two strong gamma lines. This indicates a cascade of two gamma emissions following a beta decay. The first gamma has an energy of 1173.7 keV and the second has an energy of 1332.5 keV. As seen in the decay scheme to the right, cobalt-60 occasionally decays by a high energy beta emission followed by a single gamma decay to ground state.

The presence of two different gamma energies means that there should be two different backscatter peaks and Compton edges. However, these features





Fig. 1: Decay scheme of ⁶⁰Co

were predicted to be very close to each other and are not distinguishable on the spectrum. The backscatter peak is at approximately 251 ± 20 keV, which is ~40 keV higher than expected. The Compton edge is located at about 980 ± 20 keV and is in between the two expected values.

Table

Features



Feature	Expected (keV)		
Backscatter Peak 1	210		
Backscatter Peak 2	214		
Compton Edge 1	963		
Compton Edge 2	1118		

Expected Cobalt

Spectrum

Table 2: Energy Resolution

Peak (Expected)	FWHM (keV)	Energy Resolution
1173.2 keV	187.26	15.9%
1332.5 keV	190.20	14.3%

Fig. 2: Cobalt-60 spectrum

Table 3: Cobalt Spectrum Peaks

Tuble 5. Cobait Spectium I cars					
Peak (keV)	Expected	Measured	FWHM	Centroid	# of St. Dev's Observed
Gamma 1	1173.2	1174.64 <u>+</u> 1.21	187.26	1200.2	205.8
Gamma 2	1332.5	1332.65 <u>+</u> 1.12	190.20	1344.5	250.5

The ²²Na Spectrum

In 89.8% of cases, ²²Na decays by the emission of a positron followed by the emission of a 1277 keV gamma photon to the ground state of ²²Ne. The positron annihilates with an electron, producing two gamma rays of 511 keV, the rest mass of an electron. This reflects the conversion of the masses of the positron and electron into energy by E=mc². The other 10.2% of the time, the decay begins with electron capture in the K shell, followed by the 1277 keV gamma decay. The energy resolution of the positron-electron annihilation peak is 7.6%, compared to 4.9% for the 1277 keV gamma peak.




Fig. 3: ²²Na decay scheme

Table 4: Sodium Spectrum Peaks

Peak	Expected	Measured	FWHM	Centroid	# of St. Dev's
	(keV)	(keV)			Observed
Annihilation	511	509.5 ± 0.55	37.94	499.53	604.6
Gamma	1277	1273 <u>+</u> 1.98	53.16	1277.32	264.0
Unknown	?	1820 <u>+</u> 2.1	57.8	1817	189.1




Fig. 4: Sodium-22 spectrum

The sodium-22 spectrum appeared as expected except for a high energy gamma peak of unknown origin. This peak occurred at an energy of 1820 keV, which is reasonably close to the combined energy of a sodium gamma and an annihilation gamma (1788 keV). We suspect that this high energy peak is the result of a 1277 keV gamma and 511 keV gamma entering the scintillation at the same time. This is possible given our experimental setup in which the source was place directly on the scintillation counter.

Unknown Spectra

Gamma spectroscopy can be employed to determine the constituents of an unknown radioactive sample. In this experiment, we were given two unknown sources, which we attempted to identify. We did this by comparing the unknown spectrum to the gamma energies of known sources.





Fig. 5: Unknown 1 Spectrum

Table 5: Identification of Peaks for unknown 1

Peak	Energy (keV)	# of St. Dev's Observed
A	42.8	157.0
В	875.5	391.6
C	1467.3	75.2

This unknown was determined to be a combination of ¹³⁷Cs and some kind of zinc compound. The cesium contributed Peak A, which is interpreted to be the barium x-ray peak of the spectrum, and Peak B, which is interpreted as the full energy peak. Zinc has a number of isotopes that decay with gamma energies between 1366.4 - 1560.4 keV which may be the source of Peak C. The values are not exact because adding the two spectra together cause the peaks to shift relative to one another. Also, as seen in the above experiments, the cesium calibration does not work perfectly for other spectra, which may have affected the accuracy of the calibration of Peak C.

We determined the composition of unknown 2 by adding together the spectra of various sources available in the lab until we found a reasonable match.






Fig. 6: Gamma spectrum for Unknown 2 and an example of combined spectra

We concluded that unknown 2 is composed of ¹³⁷Cs, ⁶⁰Co, ¹³³Ba, and ¹⁰⁹Cd. Peak A is the barium x-ray peak from cesium. Peaks B, C, and D are the result of combining the spectra of the four sources. They do not exactly line up with any single line of the sources since the Compton continuum from each source is in this energy region. Cobalt-60 is the source of the dual peaks E and F.

Table 6: Peak Identification for unknown 2


Peak | Energy (keV) | # of St. Dev's Observed 30.165 527.147821 A 441.722967 В 90.225 C 457.865 335.940357 843.705 775.803903 D 206.918285 \mathbf{E} 1497.09 1688.19 193.583585

Table 7: Known Energies

1	137 Cs	32 keV, 662 keV
6	⁰ Co	1173.2 keV, 1332.5 keV
1	. ³³ Ba	80.99 keV, 356.02 keV
1	09 Cd	81 keV, 356 keV

⁶⁴Cu Production and Decay

In this part of the experiment we tried to determine if gamma spectroscopy could be useful in detecting trace amounts of copper in a material of unknown composition. We exposed 16.037g of Copper to an active neutron source to create ⁶⁴Cu. This isotope has a gamma line at 1345.84 keV and a half life of 12.7 hr. The parent isotope, ⁶³Cu, has a relative abundance of 69.17%, making Cu a good candidate for this experiment.



We observed the expected peak and used the data below to calculate how many ⁶⁴Cu were produced. These "back of the envelope" calculations were done under the assumption that 20% of the emitted gammas were absorbed in the detector, which has a 12% efficiency.


Fig. 8: ⁶⁴Cu decay scheme.

Fig. 7: Irradiated copper spectrum

The signal can be approximated by a Gaussian distribution using the method of least squares.

$$G = Ae^{\frac{-(E - E_o)^2}{2\sigma^2}}$$

The above equation describes a Gaussian distribution centered at E_o with variance σ^2 . The least squares method is used to determine the constant A, as well as the error in A. We found A to be 3427 ± 175 counts. The Gaussian distribution fits well near the center of the peak, but underestimates our data toward the limbs of the peak. Assuming



the data forms a Gaussian curve, the number of standard deviations observed in the peak is 266.4.

Fig. 9: Gaussian Fit to Copper Data

Energy (keV)

◆ Copper Data ■ Gaussian Fit

We can approximate the number of ⁶⁴Cu nuclei produced while the copper was exposed to the neutron source using the data gathered with the scintillation counter and the properties of the isotope and its parent isotope, ⁶³Cu. We also assumed that 20% of the emitted gammas were absorbed in the detector, which has a 12% efficiency.

 Table 8: 64Cu Data

 Peak Energy
 1345.84 keV

 FWHM
 97.82 keV

 Net Area
 75323 Cts

 # of St. Dev
 266.4

 Sample Mass
 16.037 g

 % Abundance of 63Cu
 69.17

 Half Life of 64Cu
 12.7 hr

Sample Calculation

Nuclei =
$$16.037 \cdot 0.6917 = 11.093g = 6.68 \cdot 10^{24}$$
 amu = $5.30 \cdot 10^{22}$ nuclei $\frac{NetArea}{Time} = \frac{75323counts}{600\,\text{sec}} = 125.54counts / \text{sec}$

Half Lif & CountRate = $12.7hr \cdot 125.4counts / \text{sec} = 5739688.8counts}$
 $\frac{11479377}{0.20 \cdot 0.12} = 478307400totalcounts$

According to this calculation, there were 4.8×10^8 Copper-64 isotopes made while the sample was exposed to the neutron source. Even though this is a very small number compared to the total number of nuclei in the sample $(7.7 \times 10^{-16} \%)$ its presence was easily visible in the spectrum. The number of standard deviations observed indicates that if copper was present as 1 part in 266.4 this method could produce a detectable number of radioactive isotopes. This apparatus could therefore be used to detect trace amounts of elements with either naturally occurring radioactive isotopes, or elements that could form a radioactive isotope if it is bombarded with neutrons.

Conclusions

Gamma ray spectroscopy is an excellent method of investigating the radioactive decay of certain isotopes. Our scintillation spectroscopy setup detected the two gamma decay cascade of cobalt-60, as well as the positron-electron annihilation that characterizes the sodium-22 spectrum. Scintillation spectroscopy has greatly increased our understanding of radioactive decay since spectra have been recorded for many isotopes.

Since there is such a large base of knowledge, gamma spectroscopy is also useful in identifying unknown samples. As shown with our copper experiment, it may possibly be used to detect trace amounts of certain elements. This application should be further explored through more research.

References

Gilmore and Hemingway, <u>Practical Gamma-Ray Spectroscopy</u> Knoll, <u>Radiation Detection and Measurement</u>
Melissinos, <u>Experiments in Modern Physics</u>
Table of Isotopes, http://ie.lbl.gov/toi.html>

Appendix A

¹³⁷Cs Calibration

¹³⁷Cs has a single gamma line at 662 keV and a barium x-ray peak at 32 keV. The decay of ¹³⁷Cs starts with the emission of an electron to an excited state of ¹³⁷Ba, followed by a gamma emission in 92% of cases. The other 8% of beta emissions bring the nucleus directly to ground state. The x-ray results from electron capture in the K shell of barium followed by the emission of an x-ray, which occurs 10% of the time.

The ¹³⁷Cs spectrum shows the characteristic features of a gamma spectrum. The Compton continuum represents all the energies of 662 keV gammas that scatter in the crystal with the scattered photon escaping without depositing its full energy. The backscatter peak shows the 662 keV gammas that scatter off something outside the detector and get sent into the crystal. When lead was placed nearby the source and detector, the backscatter rate increased. The Compton continuum is a continuous range of energies because gammas can scatter at any angle from 0 to π . The highest energy of this range corresponds to the energy of the recoiling electron for a scattering angle $\theta = \pi$ in equation (3). This energy is the location of the Compton edge and the energy of the scattered photon at this scattering angle is the location of the backscattering peak.




Fig. A1: Decay scheme of ¹³⁷Cs




Fig. A2: Gamma Ray spectrum for ¹³⁷Cs with coarse gain 2, fine gain 3.

Calibration of the channels of the MCA was done for several amplifier settings and is summarized in Table 1. At every amplifier setting used, the features described above appeared at approximately the same energy once the channels were calibrated.

Table A1: Cesium Calibration

Number of Channels	CG	FG	Energy=m(channel)+b
1024	2	3	E = 3.56*C - 0.0338
2048	2	3	E = 1.82 * C - 0.775
1024	4	3	E = 1.81 * C - 2.397
1024	2	5	E = 2.395 * C - 1.536

Table 2 shows a statistical analysis of the barium x-ray and full energy peaks for an amplifier setting of CG 2 FG 3. The uncertainties were determined using the peak finding software in the MCA program and converted into keV using the calibrations calculated above.

Table A2: Cesium Spectrum Peaks

Peak	Expected	Measured	FWHM	Centroid	X ²	X ² /Degrees of
	(keV)	(keV)	(keV)	(keV)		Freedom
Barium X-	32	32.0062 <u>+</u>	27.94	24.02	3290.27	71.53
ray		0.43				
Full	662	662.1262 <u>+</u>	53.46	657.70	1095.38	52.16
Energy		0.68				

We calculated the energy resolution to be 8.2% using the full width half maximum of the full energy peak divided by the centroid of the peak.

Using equations (2) and (3) the expected energies of the Compton edge and backscatter peak were calculated. Our measured energies for these features were 7.1% and 12% different from the expected values, as shown in Table 3. However, due to the dispersed nature of these features in the spectrum, as well as calibration error, the uncertainty in the measurement was approximately 20 keV.

Table A3: Cesium Spectrum Features

(keV)	Expected	Measured	Percent Difference
Compton Edge	478	444 <u>+</u> 20	7.1%
Backscatter Peak	184	206 <u>+</u> 20	12%

Appendix B

Notes on Sources

There is usually more than one decay path for a radioactive source, which have different probabilities for occurrence. The half lives of intermediate states may also have different half lives. Due to these differences, the multiple gamma peaks of certain radioactive nuclides may vary in relative strength as time passes from the preparation of the source. The following table lists the characteristics of the known sealed sources used.

Table B1: Source Data

Source	Activity (μCi)	Half Life	Prep Date
¹³⁷ Cs	1.0	30.2 y	June 2005
⁶⁰ Co	1.0	6.27 y	Sept 2000
²² Na	1.0	2.6	June 2005

Appendix C

Energy Resolution

The energy resolution of the scintillator is described by the fractional width of a peak in a spectrum. The width of a peak is related to the number of counts in a peak, which is a function of the number of visible photons that are generated and detected (N).

$$\Delta E/E = \Delta N/N = 1/\sqrt{N} \alpha 1/\sqrt{E}$$
 (5)

Equation (5) shows the relation between fractional width and the probability of photon production in the scintillator. For 137 Cs, the number of visible photons produced and detected (N) is estimated to be 148. The width of a peak should therefore scale as a function of $1/\sqrt{N}$. From this we can calculated the expected FWHM of the Cesium gamma peak using equation (6).

$$\Delta E = E/\sqrt{N}$$
 (6)

Using E=662 and N=148, the expected FWHM for ¹³⁷Cs is 54.4 keV. The measured value is 53.46 keV, which differs from the expected by 1.7%.

The fractional width should also scale linearly with $1/\sqrt{E}$. The plot below shows that more data for intermediate energies is necessary to verify this relation.




Fig. C1: Fractional width vs. $1/\sqrt{E}$ with linear trend lines. Purple line shows trend for all gamma energies, blue line shows gamma energies for Cesium full energy, Sodium and Copper peaks.